

Instruments & Implants (I&I)

Precision machining of complex orthopedic instruments, complete instrument sets, and implant components and assemblies

Design | Develop | Manufacture | Package | Deliver

Why Choose **Avalign** I&I?

- Expertise in implant components and assemblies, complex instruments, and full kitting across all segments of orthopedics
- Latest CNC technology across swiss lathe, mill, mill-turn, grinding, wire EDM, and gun drill
- Blast, polish, finish, clean and passivate, laser etch ... all in one location
- Advanced technology, including CT Metrology and our proprietary Flex Technology
- Significant new capacity installed and available to support your next project



Your Trusted CDMO Partner

Avalign is one of the only CDMOs that can provide every aspect of an orthopedic implant system, including: Implants and Coatings | Fixation Devices | Complex Instruments | Cutting Instruments | German Specialty Instruments | Cases & Trays | Fully Kitted Implant Banks and Instrument Sets

98%

On Time
Delivery
Commitment

98%

Customer Lot
Acceptance
Commitment

25+

Years
Experience

12-16

Week Lead
Time

Contact Us



(855) 282-5446



www.avalign.com